		Revisions	SO 0191-5
Rev from	Rev to	ECO #	
0602	0404	2771-04	
2. Blank (Numbe Page S	er of Pages: 16	ength: 8.5" Width: 11" Number of Sheets: 8 " Width 5.5" Final Folded S	Size: 8.5" x 5.5"
	#1	#2 #3	#4
	HEADER #5	HEADER #6 #7	
5. Ink Col No. of 6. Graphi	lors: Printed tv Colors: 1 ics are approve	otrol Number 8809911 for Material Inf vo sides X Yes No PMS# 513 Purple ed by Becton, Dickinson and Company. ng the most current approved revision	Supplier has the
abel Design roofer hecked By	Date Date Date	COMPANY CONFIDENTAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION Category and Description	Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA Sheet: 1 of 13
art Number: 88	309911JAA	Package Insert, BBL Crystal GP ID Kit	Scale: N/A

Solution Systems Gram-Positive ID Kit

8809911JAA 2004/04 CLIA COMPLEXITY: HIGH CDC IDENTIFIER CODES ANALYTE: 0412 TEST SYSTEM: 07919

> U.S. Pat. 5,182,082 U.S. Pat. 5,338,666

See symbol glossary at end of insert.

INTENDED USE

The **BBLCrystal™** Gram-Positive (GP) Identification (ID) system is a miniaturized identification method employing modified conventional, fluorogenic and chromogenic substrates. It is intended for the identification of aerobic gram-positive bacteria.^{1,2,13,16}

SUMMARY AND EXPLANATION

Micromethods for the biochemical identification of microorganisms were reported as early as 1918.³ Several publications reported on the use of the reagent-impregnated paper discs and micro-tube methods for differentiating enteric bacteria.^{3,4,7,17,19} The interest in miniaturized identification systems led to the introduction of several commercial systems in the late 1960s, and they provided advantages in requiring little storage space, extended shelf life, standardized quality control and ease of use.

In general, many of the tests used in the **BBLCrystal** ID Systems are modifications of classical methods. These include tests for fermentation, oxidation, degradation and hydrolysis of various substrates. In addition, there are chromogen and fluorogen linked substrates, as in the **BBLCrystal** GP ID panel, to detect enzymes that microbes use to metabolize various substrates.^{5,7,8,9,11,12,14,15}

The **BBLCrystal** GP ID kit is comprised of (i) **BBLCrystal** GP ID panel lids, (ii) **BBLCrystal** bases and (iii) **BBLCrystal** ANR, GP, RGP, N/H ID Inoculum Fluid (IF) tubes. The lid contains 29 dehydrated substrates and a fluorescence control on tips of plastic prongs. The base has 30 reaction wells. Test inoculum is prepared with the inoculum fluid and is used to fill all 30 wells in the base. When the lid is aligned with the base and snapped in place, the test inoculum rehydrates the dried substrates and initiates test reactions.

Following an incubation period, the wells are examined for color changes or presence of fluorescence that result from metabolic activities of the microorganisms. The resulting pattern of the 29 reactions is converted into a ten-digit profile number that is used as the basis for identification.¹⁸ Biochemical and enzymatic reaction patterns for the 29 **BBL Crystal** GP ID substrates for a wide variety of microorganisms are stored in the **BBL Crystal** GP ID data base. Identification is derived from a comparative analysis of the reaction pattern of the test isolate to those held in the database. A complete list of taxa that comprises the current database is provided in **Table 1** (see pg. 7).

PRINCIPLES OF THE PROCEDURE

The **BBLCrystal** GP ID panels contain 29 dried biochemical and enzymatic substrates. A bacterial suspension in the inoculum fluid is used for rehydration of the substrates. The tests used in the system are based on microbial utilization and degradation of specific substrates detected by various indicator systems. Enzymatic hydrolysis of fluorogenic substrates containing coumarin derivatives of 4-methylumbelliferone (4MU) or 7-amino-4-methylcoumarin (7-AMC), results in increased fluorescence that is easily detected visually with a UV light source.^{11,12,14,15} Chromogenic substrates upon hydrolysis produce color changes that can be detected visually. In addition, there are tests that detect the ability of an organism to hydrolyze, degrade, reduce or otherwise utilize a substrate in the **BBLCrystal** ID Systems.

Reactions employed by various substrates and a brief explanation of the principles employed in the system are described in **Table 2** (see pg. 8). Panel location in referred tables indicates the row and column where the well is located (example: 1J refers to Row 1 in column J).

REAGENTS

The **BBLCrystal** GP ID panel contains 29 enzymatic and biochemical substrates. Refer to **Table 3** (see pg. 9) for a list of active ingredients.

Warnings and Precautions:

For in vitro Diagnostic Use.

After review by the U.S. Centers for Disease Control and Prevention (CDC), and the Food and Drug Administration (FDA) under CLIA '88, this product has been identified as high complexity. The CDC Analyte Identifier Code is 0412; the CDC Test System Identifier Code is 07919.

After use, all infectious materials including plates, cotton swabs, inoculum fluid tubes, and panels must be autoclaved prior to disposal or incineration.

STORAGE AND HANDLING/SHELF LIFE

Lids: BBL Crystal GP lids are individually packaged and must be stored unopened in a refrigerator at 2–8°C. DO NOT FREEZE. Visually inspect the package for holes or cracks in the foil package. Do not use if the packaging appears to be damaged. Lids in the original packaging, if stored as recommended, will retain expected reactivity until the date of expiration.

Bases: Bases are packaged in two sets of ten, in **BBL Crystal** incubation trays. The bases are stacked facing down to minimize air contamination. Store in a dust-free environment at 2–30°C, until ready to use. Store unused bases in the tray, in plastic bag. Empty trays should be used to incubate inoculated panels.

Inoculum Fluid: BBL Crystal ANR, GP, RGP, N/H ID Inoculum Fluid (IF) is packaged in two sets of ten tubes. Visually inspect the tubes for cracks, leaks, etc. Do not use if there appears to be a leak, tube or cap damage or visual evidence of contamination (i.e., haziness, turbidity). Store tubes at 2–25°C. Expiration dating is shown on the tube label. Only ANR, GP, RGP, N/H Inoculum Fluid should be used with **BBL Crystal** GP ID panels.

On receipt, store the **BBL Crystal** GP ID kit at 2–8°C. Once opened, only the lids need to be stored at 2–8°C. The remaining components of the kit may be stored at 2–25°C. If the kit or any of the components are stored refrigerated, each should be brought to room temperature prior to use.

SPECIMEN COLLECTION AND PROCESSING

BBL Crystal ID Systems are not for use directly with clinical specimens. Use isolates from media such as Trypticase™ Soy Agar with 5% Sheep Blood (TSA II) or Columbia Agar with 5% Sheep Blood (Columbia Blood Agar). Use of selective media such as Phenylethyl Alcohol Agar with 5% Sheep Blood (PEA) or Columbia CNA Agar with 5% Sheep Blood (CNA) is also acceptable. Media containing esculin should not be used. The test isolate must be a pure culture, no more than 18–24 h old for most genera; for some slow growing organisms up to 48 h may be acceptable. When swabs are utilized, only cotton-tipped applicators should be used to prepare the inoculum suspensions. Some polyester swabs may cause problems with inoculation of the panels. (See "Limitations of the Procedure").

The incubator used should be humidified to prevent evaporation of fluid from the wells during incubation. The recommended humidity level is 40-60%. The usefulness of **BBLCrystal** ID Systems or any other diagnostic procedure performed on clinical specimens is directly influenced by the quality of the specimens themselves. It is strongly recommended that laboratories employ methods discussed in the *Manual of Clinical Microbiology* for specimen collection, transport and inoculation onto primary isolation media.^{1,16}

TEST PROCEDURE

Materials Provided: BBL Crystal GP ID Kit -

20 BBL Crystal GP ID Panel Lids,

20 BBL Crystal Bases,

20 **BBL Crystal** ANR, GP, RGP, N/H ID IF Tubes. Each tube has approximately 2.3 ± 0.15 mL of Inoculum Fluid containing: KCl 7.5 g, CaCl₂ 0.5 g, Tricine N-[2-Hydroxy-1, 1-bis (hydroxymethyl)methyl] glycine 0.895 g, purified water to 1000 mL,


2 incubation trays,

1 BBL Crystal GP ID Color Reaction Chart and Results Pad.

Materials Required But Not Provided: Sterile cotton swabs (*do not use polyester swabs*), incubator (35 – 37°C) non-CO₂ (40–60% humidity), McFarland No. 0.5 standard, **BBL Crystal** Panel Viewer, **BBL Crystal** ID System Electronic Codebook or **BBL Crystal** Gram-Positive Manual Codebook, and appropriate culture media. Also required are the necessary equipment and labware used for preparation, storage and handling of clinical specimens.


Test Procedure: BBL Crystal GP ID System requires a Gram stain.

- Remove lids from pouch. Discard desiccant. Once removed from the pouch, covered lids should be used within 1 h. Do not use the panel if there is no desiccant in the pouch.
- Take an inoculum fluid tube and label with patient's specimen number. Using aseptic technique, pick colonies of the same morphology with the tip of a sterile cotton swab (do not use a polyester swab) or a wooden applicator stick from one of the recommended media (see section under "Specimen Collection and Processing").
- Suspend colonies in a tube of BBL Crystal ANR, GP, RGP, N/H ID Inoculum Fluid.




- Recap tube and vortex for approximately 10–15 s. The turbidity should be equivalent to a McFarland No. 0.5 standard. If the inoculum suspension concentration is in excess of the recommended McFarland standard, one of the following steps is recommended:
 - a. Use a fresh tube of inoculum fluid to prepare a new inoculum suspension equivalent to a McFarland No. 0.5 standard.
 - b. If additional colonies are unavailable for preparation of a new inoculum suspension, using aseptic techniques, dilute the inoculum by adding the minimum required volume (not to exceed 1.0 mL) of 0.85% sterile saline or inoculum fluid to bring down the turbidity equivalent to a McFarland No. 0.5 standard. Remove the excess amount added to the tube with a sterile pipet so that the final volume of inoculum fluid is approximately equivalent to that of the original volume in the tube (2.3 ± 0.15 mL). Failure to make this adjustment in volume will result in spilling of the inoculum suspension over the black portion of the base rendering the panel unusable.

- 5. Take a base, and mark the patient's specimen number on the side wall.
- 6. Pour entire contents of the inoculum fluid tube into target area of the base.



- 7. Hold base in both hands and roll inoculum gently along the tracks until all of the wells are filled. Roll *back* any excess fluid to the target area and place the base on a bench top.
- 8. Align the lid so that the labeled end of the lid is on top of the target area of the base.



 Push down until a slight resistance is felt. Place thumb on edge of lid towards middle of panel on each side and push downwards simultaneously until the lid snaps into place (listen for two "clicks").

Purity Plate: Using a sterile loop, recover a small drop from the inoculum fluid tube either before or after inoculating the base and inoculate an agar slant or plate (any appropriate medium) for purity check. Discard inoculum fluid tube and cap in a biohazard disposal container. Incubate the slant or plate for 24–48 h at 35–37°C under appropriate conditions. The purity plate or slant may also be used for any supplementary tests or serology, if required.

Incubation: Place inoculated panels in incubation trays. Ten panels can fit in one tray (5 rows of 2 panels). All panels should be incubated **face down** (larger windows facing up; label facing down) in a non-CO₂ incubator with 40–60% **humidity**. Trays should not be stacked more than two high during incubation. The incubation time for panels is 18–24 h at 35–37°C. If panels are incubated for 24 h, they should be read within 30 min after removing from incubator.

Reading: After the recommended period of incubation, remove the panels from the incubator. All panels should be read **face down**


panels from the incubator. All panels should be read face down (larger windows up; label facing down) using the **BBLCrystal** Panel Viewer. Refer to the color reaction chart and/or **Table 3** (see pg. 9) for an interpretation of the reactions. Use the results pad to record reactions. Alternatively, the **BBLCrystal** AutoReader may be used to read the panels.

- a. Read columns E thru J first, using the regular (white) light source.
- b. Read columns A thru D (fluorescent substrates) using the UV light source in the panel viewer. A fluorescent substrate well is considered positive *only* if the intensity of the fluorescence observed in the well is *greater* than the Negative Control well (4A).

Calculation of BBL Crystal Profile Number: Each test result (except 4A, which is used as a fluorescence negative control) scored positive is given a value of 4, 2, or 1, corresponding to the row where the test is located. A value







of 0 (zero) is given to any negative result. The values resulting from each positive reaction in each column are then added together. A 10-digit number is generated; this is the profile number.

Example:	Α	В	С	D	E	F	G	н	I	L
4	*	+	-	-	+	+	+	-	+	-
2	-	+	+	+	-	+	-	+	+	-
1	+	-	+	-	+	-	-	+	+	-
Profile	1	6	3	2	5	6	4	3	7	0

*(4A) = fluorescent negative control

The resulting profile number and cell morphology, if known, should be entered on a PC in which the **BBL Crystal** Mind Software has been installed to obtain the identification. If using the **BBL Crystal** AutoReader, organisms are automatically identified by the PC. A manual codebook is also available. If a PC is not available contact BD Technical Services for assistance with the identification.

User Quality Control: Quality control testing is recommended for each lot of panels as follows -

- 1. Inoculate a panel with Streptococcus pyogenes ATCC™ 19615 per recommended procedure (refer to 'Test Procedure")
- Incubate panel for 18-20 h at 35-37°C.
- Read panel with the panel viewer and color reaction chart; record reactions using the results pad. 3. Alternatively, read the panel on the BBL Crystal AutoReader.
- Compare recorded reactions with those listed in **Table 4** (see *pg. 10*). If discrepant results are obtained, confirm purity of quality control strain before contacting BD Technical Services. 4.

Expected test results for additional quality control test strains are listed in Table 5 (see pg. 10). Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent NCCLS guidance and CLIA regulations for appropriate Quality Control practices.

LIMITATIONS OF THE PROCEDURE

The BBL Crystal GP ID System is designed for the taxa provided. Taxa other than those listed in Table 1 are not intended for use in this system.

The BBL Crystal GP ID System database includes some species that are rarely isolated from human clinical specimens and were not encountered in the clinical studies of this product. It also includes some species that were encountered less than 10 times in the clinical studies. Refer to **Table 1** (see pg. 7) for a breakdown of the number of strains per species tested in clinical trials. The laboratorian should determine if additional testing is required to confirm identity of those species for which performance has not been established (i.e., those species where less than 10 isolates were evaluated in the clinical trials for this product).

The BBL Crystal GP ID database was developed with BBL™ brand media. Reactivity of some substrates in miniaturized identification systems may be dependent upon the source media used in inoculum preparations. We recommend the use of the following media for use with the **BBL Crystal** GP ID System: TSA II and Columbia Blood Agar. Use of selective media, such as PEA or CNA is also acceptable. Media containing esculin should not be used. BBL Crystal Identification Systems use a modified microenvironment; therefore, expected values for its individual tests may differ from information previously established with conventional test reactions. The accuracy of the

BBL Crystal GP ID System is based on statistical use of specially designed tests and an exclusive database. While BBL Crystal GP ID System aids in microbial differentiation, it should be recognized that minor variations may exist in strains within species. Use of panels and interpretation of results require a competent microbiologist. The final identification of the isolate should take into consideration the source of the specimen, aerotolerance, cell morphology, colonial characteristics on various media as well as metabolic end products as determined by gas-liquid chromatography, when warranted.

Only cotton-tipped applicator swabs should be used to prepare the inoculum suspension as some polyester swabs may cause the inoculum fluid to become viscous. This may result in insufficient inoculum fluid to fill the wells. Covered lids once removed from the sealed pouches must be used within 1 h to ensure adequate performance. The incubator where panels are placed should be humidified to prevent evaporation of inoculum fluid from the

wells during incubation. The recommended humidity level is 40-60%. The panels, after inoculation, should only be incubated face down (larger windows facing up; label facing down)

to maximize the effectiveness of substrates.

If the BBL Crystal test profile yields a "No identification" result and culture purity has been confirmed, then it is likely that (i) the test isolate is producing atypical BBL Crystal reactions (which may also be caused by procedural errors), (ii) the test species is not part of the intended taxa or (iii) the system is unable to identify the test isolate with the required level of confidence. Conventional test methods are recommended when user error has been ruled out

EXPECTED VALUES

The expected substrate reactions for the species of organisms most frequently encountered in the clinical study of **BBL Crystal** GP ID System are shown in **Table 6** (see *pg. 11*). The provided percentages were generated from reactions given by the organisms used in generating the database. **Table 1** (see *pg. 7*) shows all the taxa tested during database generation.

PERFORMANCE CHARACTERISTICS

Reproducibility: In an external study involving four clinical laboratories, (total of four evaluations), the reproducibility of **BBLCrystal** GP ID substrates' (29) reactions was studied by replicate testing. The reproducibility of the individual substrate reactions ranged from 79.2%–100%. The overall reproducibility of **BBLCrystal** GP ID panel was determined to be 96.7%.²⁰

Accuracy of Identification: The performance of BBL Crystal GP ID System was compared to currently available commercial systems using <u>clinical isolates and stock cultures</u>. A total of four studies were conducted in four independent laboratories. Fresh, routine isolates arriving in the clinical laboratory, as well as previously identified isolates of the clinical trial sites' choice, were utilized to establish performance characteristics. Out of 735 total isolates tested from the four studies using BBL Crystal GP Identification System, 623 (84.8%) were correctly identified without the use of supplement tests, and 668 (90.9%) were correctly identified when supplemental tests were included. A total of 56 (7.6%) isolates were incorrectly identified, and a message of "No Identification" was obtained for 11 (1.5%) isolates.²⁰ Table 7 (see pg. 12) shows the accuracy of

identification for the species most frequently encountered (i.e., 10 or more isolates) in the clinical trial as well as for the remaining group of species where less than 10 isolates were tested.

Cat. No.

Description

AVAILABILITY

Cat. No. Description

	2 coulput		Description
245240	BBL Crystal™ Gram-Positive ID Kit,	245300	BBLCrystal™ AutoReader
	containing 20 each: BBL Crystal GP ID Panel Lids, BBL Crystal Bases and	221165	BBL™ Columbia Agar with 5% Sheep Blood, pkg. of 20.
	BBL Crystal ANR, GP, RGP, N/H ID Inoculum Fluid.	221263	BBL™ Columbia Agar with 5% Sheep Blood, ctn. of 100.
245038	BBL Crystal™ ANR, GP, RGP, N/H ID Inoculum Fluid, ctn. of 10.	221352	BBL™ Columbia CNA Agar with 5% Sheep Blood, pkg. of 20.
245031	BBL Crystal™ Panel Viewer, Domestic model, 110 V, 60 Hz.	221353	BBL™ Columbia CNA Agar with 5% Sheep Blood, ctn. of 100.
245032	BBL Crystal™ Panel Viewer, European model, 220 V, 50 Hz.	221179	BBL™ Phenylethyl Alcohol Agar with 5% Sheep Blood, pkg. of 20.
245033	BBL Crystal™ Panel Viewer, Japanese model, 100 V, 50/60 Hz.	221277	BBL™ Phenylethyl Alcohol Agar with 5% Sheep Blood, ctn. of 100.
245034	BBLCrystal™ Panel Viewer, Longwave UV Tube.	221239	BBL™ Trypticase™ Soy Agar with 5% Sheep Blood (TSA II), pkg. of 20.
245036	BBLCrystal™ Panel Viewer, White Light Tube.	221261	BBL™ Trypticase™ Soy Agar with 5% Sheep Blood (TSA II), ctn. of 100.
245037	BBL Crystal™ Identification Systems Gram-Positive Manual Codebook.	212539	BBL™ Gram Stain Kit, pkg. of 4 x 250 mL bottles.
441010	BBL Crystal [™] Mind Software		bottles.

REFERENCES

- Balows, A., W.J. Hausler, Jr., K.L. Herrmann, H.D. Isenberg, and H.J. Shadomy (ed.). 1991. Manual of clinical microbiology, 5th ed. American Society for Microbiology, Washington, D.C.
- Baron, E.J., L.R. Peterson, and S.M. Finegold. 1994. Bailey and Scott's diagnostic microbiology, 9th ed. Mosby-Year Book, Inc., St. Louis.
- Bronfenbrenner, J., and M.J. Schlesinger. 1918. A rapid method for the identification of bacteria fermenting carbohydrates. Am. J. Public Health. 8:922-923.
- Cowan, S.T., and K.J. Steel. 1974. Manual for the identification of medical bacteria. 2nd ed. Cambridge University Press, Cambridge.
- Edberg, S.C., and C.M. Kontnick. 1986. Comparison of β-glucuronidase-based substrate systems for identification of *Escherichia coli*. J. Clin. Microbiol. 24:368-371.
- Ferguson, W.W., and A.E. Hook. 1943. Urease activity of *Proteus* and *Salmonella* organisms. J. Lab. Clin. Med. 28:1715-1720.
- 7. Hartman, P.A. 1968. Miniaturized microbiological methods. Academic Press, New York.
- Kampfer, P., O. Rauhoff, and W. Dott. 1991. Glycosidase profiles of members of the family Enterobacteriaceae. J. Clin. Microbiol. 29:2877-2879.
- 9. Killian, M., and P. Bulow. 1976. Rapid diagnosis of *Enterobacteriaceae* 1: detection of bacterial glycosidases. Acta Pathol. Microbiol. Scand. Sect. B. 84:245-251.
- 10. MacFaddin, J.F. 1980. Biochemical tests for identification of medical bacteria, 2nd ed. Williams & Wilkins, Baltimore.
- 11. Maddocks, J.L., and M. Greenan. 1975. Rapid method for identifying bacterial enzymes. J. Clin. Pathol. 28:686-687. 12. Manafi, M., W. Kneifel, and S. Bascomb. 1991. Fluorogenic and chromogenic substrates used in bacterial
- diagnostics. Microbiol. Rev. 55:335-348.
- Mandell, G.L., R.G. Douglas, Jr. and J.E. Bennett. 1990. Principles and practice of infectious diseases, 3rd ed. Churchill Livingstone Inc., New York.
- 14. Mangels, J., I. Edvalson, and M. Cox. 1993. Rapid Identification of *Bacteroides fragilis* group organisms with the use of 4-methylumbelliferone derivative substrates. Clin. Infect. Dis. *16*(54):5319-5321.

- Moncla, B.J., P. Braham, L.K. Rabe, and S.L. Hiller. 1991. Rapid presumptive identification of black-pigmented gram-negative anaerobic bacteria by using 4-methylumbelliferone derivatives. J. Clin. Microbiol. 29:1955-1958.
- Murray, P.R., E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Yolken (ed.). 1995. Manual of clinical microbiology, 6th ed. American Society for Microbiology, Washington, D.C.
- Sanders, A.C., J.E. Faber, and T.M. Cook. 1957. A rapid method for the characterization of enteric pathogen using paper discs. Appl. Microbiol. 5:36-40.
- 18. Sneath, P.H.A. 1957. The application of computers to taxonomy. J. Gen. Microbiol. 17:201-221.
- 19. Soto, O.B. 1949. Fermentation reactions with dried paper discs containing carbohydrate and indicator. Puerto Rican J. Public Health. Trop. Med. 25:96-100.
- 20. Data on file at BD Diagnostics.

Key: Table 1

KEY:

* = These taxa have fewer than 10 unique BBL Crystal profiles in the current database.

- ("x") = Number of isolates (i.e., "x") encountered in the clinical trial. If no number in parenthesis is shown after an organism name or group description, these species were not encountered in the clinical trial.
- Note #1: There were 14 additional isolates encountered in the clinical trial that are not shown above. Five (5) (i.e., 4 Staphylococcus species and 1 Enterococcus were identified only to the genus level by the reference system against which BBLCrystal GP was compared, although BBLCrystal GP identified these organisms to the species level. Nine (9) were identified by the reference system, but were not included in the BBLCrystal GP database taxa.

Note #2: The organisms shown in bold face type were encountered 10 or more times in the clinical study for this product.

Note #3: The organisms not shown in bold face type are either species which are rarely isolated from human clinical specimens or species that were infrequently (less than 10) encountered in the clinical study for this product. The laboratorian should determine if additional testing is required to confirm their identity.

Table 1

Taxa in BBLCrystal™ GP ID System

Actinomyces pyogenes Aerococcus species (includes A. urinae and A. viridans) Aerococcus urinae Aerococcus viridans Alloiococcus otitidis* Arcanobacterium haemolyticum*(2) Bacillus brevis (1) Bacillus cereus (2) Bacillus circulans Bacillus coagulans Bacillus licheniformis (1) Bacillus megaterium Bacillus pumilus Bacillus species (includes B. brevis, B. circulans, B. coagulans, B. licheniformis, B. megaterium, B. pumilus and B. sphaericus, P. alvei, P. macerans) (9) Bacillus sphaericus Bacillus subtilis (1) Corvnebacterium aquaticum Corynebacterium bovis Corynebacterium diphtheriae (includes C. diphtheriae subsp gravis, C. diphtheriae subsp mitis and C diphtheriae subsp intermedius) Corynebacterium genitalium Corynebacterium jeikeium (7) Corynebacterium kutscheri Corynebacterium propinquum (1) Corvnebacterium pseudodiphtheriticum (2) Convnebacterium pseudogenitalium Corvnebacterium pseudotuberculosis (2) Corynebacterium renale group Corvnebacterium species (includes C. aquaticum, C. bovis, C. kutscheri, C. propinguum, C. pseudodiphtheriticum, C. pseudotuberculosis C. renale group, C. striatum and C. ulcerans) (29) Corynebacterium striatum (6) Corynebacterium ulcerans Enterococcus avium (3) Enterococcus casseliflavus/gallinarum (14)

Enterococcus durans (2) Enterococcus faecalis (78) Enterococcus faecium (33) Enterococcus hirae

Enterococcus raffinosus (3) Enterococcus solitarius Erysipelothrix rhusiopathiae Gardnerella vaginalis Gemella haemolysans Gemella morbillorum Gemella species (includes G. haemolysans and G. morbillorum) Globicatella sanguis (3) Helcococcus kunzii Lactococcus garvieae Lactococcus lactis subsp cremoris Lactococcus lactis subsp hordniae Lactococcus lactis subsp lactis Lactococcus raffinolactis Lactococcus species (includes L. lactis subsp cremoris, L. lactis subsp hordniae, L. lactis subsp lactis and L. raffinolactis) Leuconostoc citreum Leuconostoc lactis (1) Leuconostoc mesenteroides subsp mesenteroides Leuconostoc pseudomesenteroides Leuconostoc species (includes L. citreum, L. lactis, L. mesenteroides subsp mesenteroides and L. pseudomesenteroides) Listeria grayi * Listeria ivanovii subsp ivanovii Listeria monocytogenes (3) Listeria murrayi Micrococcus kristinae Micrococcus luteus Micrococcus lylae Micrococcus roseus Micrococcus sedentarius Micrococcus species (includes M. kristinae, M. luteus, M. lylae, M. roseus and M. sedentarius) (10) Oerskovia species (includes O. turbata and O. xanthineolytica)

Paenibacillus alvei

Pediococcus species (includes P. damnosus, P. parvulus and P. pentosaceus) Rhodococcus equi Rothia dentocariosa* (1) Staphylococcus aureus (88) Staphylococcus auricularis (2) Staphylococcus capitis (includes *S. capitis* subsp capitis and *S. capitis* subsp ureolyticus) (13) Staphylococcus caprae Staphylococcus carnosus Staphylococcus cohnii (includes S. cohnii subsp cohnii and S. cohnii subsp urealyticum) (1) Staphylococcus cohnii subsp cohnii Staphylococcus cohnii subsp urealyticum Staphylococcus epidermidis (88) Staphylococcus equorum Staphylococcus felis Staphylococcus gallinarum Staphylococcus haemolyticus (23) Staphylococcus hominis (17) Staphylococcus intermedius Staphylococcus kloosii (2) Staphylococcus lentus Staphylococcus lugdunensis (3) Staphylococcus pasteuri *(1) Staphylococcus saccharolyticus (6) Staphylococcus saprophyticus Staphylococcus schleiferi (includes S. schleiferi subsp coagulans and S. schleiferi subsp schleiferi) Staphylococcus sciuri Staphylococcus simulans (3)

Paenibacillus macerans

Pediococcus damnosus

Pediococcus parvulus

Pediococcus pentosaceus

Staphylococcus vitulus Staphylococcus warneri (6) Staphylococcus xylosus (1) Stomatococcus mucilaginosus (6)

Streptococcus acidominimus

Streptococcus agalactiae (54)

Streptococcus anginosus (1) Streptococcus bovis (includes S. bovis I and S. bovis II) (10) Streptococcus constellatus (1) Streptococcus cricetus * Streptococcus crista Streptococcus equi (includes S. equi subsp equi and S. equi subsp zooepidemicus) (1) Streptococcus equi subsp equi (2) Streptococcus equi subsp zooepidemicus Streptococcus equinus Streptococcus gordonii Streptococcus Group C/G (11) Streptococcus intermedius Streptococcus milleri group (includes S. anginosus, S constellatus and S. intermedius) (20) Streptococcus mitis (4) Streptococcus mitis group (includes S. mitis and . S. oralis) (20) Streptococcus mutans Streptococcus mutans group (includes S. cricetus. S. mutans and S. sobrinus) (2) Streptococcus oralis Streptococcus parasanguis (1) Streptococcus pneumoniae (54) Streptococcus porcinus Streptococcus pyogenes (50) Streptococcus salivarius (3) Streptococcus salivarius group (includes S. salivarius and S. vestibularis) (4)

Streptococcus sanguis (2) Streptococcus sanguis group (includes S. crista, S. gordonii, S. parasanguis and S. sanguis) Streptococcus sobrinus Streptococcus uberi Streptococcus vestibularis Turicella otitidis *

Panel Location	Test Feature	Code	Principle (Reference)
4A	Fluorescent negative control	FCT	Control to standardize fluorescent substrate results.
2A	4MU-β-D-glucoside	FGC	-
1A	L-valine-AMC	FVA	_
4B	L-phenylalanine-AMC	FPH	_
2B	4MU-α-D-glucoside	FGS	-
1B	L-pyroglutamic acid-AMC	FPY	Enzymatic hydrolysis of the amide or glycosidic
4C	L-tryptophan-AMC	FTR	bond results in the release of a fluorescent coumarin derivative. ^{5,8,11,12,14,15}
2C	L-arginine-AMC	FAR	-
1C	4MU-N-acetyl-β-D-glucosaminide	FGA	-
4D	4MU-phosphate	FHO	_
2D	4MU-β-D-glucuronide	FGN	
1D	L-isoleucine-AMC	FIS	-
4E	Trehalose	TRE	_
2E	Lactose	LAC	
1E	Methyl- α & β -glucoside	MAB	_
4F	Sucrose	SUC	-
2F	Mannitol	MNT	Utilization of carbohydrate results in lower pH and
1F	Maltotriose	MTT	change in indicator (Phenol red). ^{1,2,3,4,7,16}
4G	Arabinose	ARA	-
2G	Glycerol	GLR	-
1G	Fructose	FRU	-
4H	p-nitrophenyl-β-D-glucoside	BGL	Enzymatic hydrolysis of the colorless aryl substituted glycoside releases yellow p-nitrophenol. ^{5,9,12}
2H	p-nitrophenyl-β-D-cellobioside	PCE	
1H	Proline & Leucine-p-nitroanilide	PLN	Enzymatic hydrolysis of the colorless amide substrate releases yellow p-nitroaniline. ^{5,9,12}
41	p-nitrophenyl-phosphate	PHO	
21	p-nitrophenyl- α -D-maltoside	PAM	Enzymatic hydrolysis of the colorless aryl substituted glycoside releases yellow p-nitrophenol. ^{5,9,12}
11	o-nitrophenyl-β-D-galactoside (ONPG) & p-nitrophenyl-α-D-galactoside	PGO	-
4J	Urea	URE	Hydrolysis of urea and the resulting ammonia change the pH indicator color (Bromthymol blue). ^{2,6,10}
2J	Esculin	ESC	Hydrolysis of esculin results in a black precipitat in the presence of ferric ion. ¹⁰
1J	Arginine	ARG	Utilization of arginine results in pH rise and change in the color of the indicator (Bromcreso purple). ²

Table 2 Principles of Tests Employed in the BBLCrystal™ GP ID System

Table 3	
Reagents used in the BBL Crystal™ GP ID System	

Panel Location	Substrate	Code	Pos.	Neg.	Active Ingredients	Approx. Amt (g/L)
4A	Fluorescent negative control	FCT	n/a	n/a	Fluorescent coumarin derivative	≤¢
2A	4MU-β-D-glucoside	FGC	blue fluorescence >FCT well	blue fluorescence ¢FCT well	4MU-β-D-glucoside	≤¢
1A	L-valine-AMC	FVA	blue fluorescence >FCT well	blue fluorescence ¢FCT well	L-valine-AMC	ţ1
4B	L-phenylalanine-AMC	FPH	blue fluorescence >FCT well	blue fluorescence ¢FCT well	L-phenylalanine-AMC	ţ1
2B	4MU-α-D-glucoside	FGS	blue fluorescence >FCT well	blue fluorescence ≤ ĘCT well	4MU-α-D-glucoside	≤¢
1B	L-pyroglutamic acid-AMC	FPY	blue fluorescence >FCT well	blue fluorescence ≤ ĘCT well	L-pyroglutamic acid-AMC	≤¢
4C	L-tryptophan-AMC	FTR	blue fluorescence >FCT well	blue fluorescence ≤ ¢ CT well	L-tryptophan-AMC	¢1
2C	L-arginine-AMC	FAR	blue fluorescence >FCT well	blue fluorescence ≤ ¢ CT well	L-arginine-AMC	≤¢
1C	4MU-N-acetyl-β-D- glucosaminide	FGA	blue fluorescence >FCT well	blue fluorescence ≤ ¢ CT well	4MU-N-acetyl-β-D- glucosaminide	≤¢
4D	4MU-phosphate	FHO	blue fluorescence >FCT well	blue fluorescence ≤ ¢ CT well	4MU-phosphate	≤¢
2D	4MU-β-D-glucuronide	FGN	blue fluorescence >FCT well	blue fluorescence ≤€CT well	4MU-β-D-glucuronide	≤¢
1D	L-isoleucine-AMC	FIS	blue fluorescence >FCT well	blue fluorescence ≤ ¢ CT well	L-isoleucine-AMC	≤¢
4E	Trehalose	TRE	Gold/Yellow	Orange/Red	Trehalose	≤ ≩ 00
2E	Lactose	LAC	Gold/Yellow	Orange/Red	Lactose	≤ ≩ 00
1E	Methyl-α & β-glucoside	MAB	Gold/Yellow	Orange/Red	Methyl-α & β-glucoside	≤ ≩ 00
4F	Sucrose	SUC	Gold/Yellow	Orange/Red	Sucrose	¢300
2F	Mannitol	MNT	Gold/Yellow	Orange/Red	Mannitol	¢300
1F	Maltotriose	MTT	Gold/Yellow	Orange/Red	Maltotriose	≤≩00
4G	Arabinose	ARA	Gold/Yellow	Orange/Red	Arabinose	¢300
2G	Glycerol	GLR	Gold/Yellow	Orange/Red	Glycerol	≤ ₿ 00
1G	Fructose	FRU	Gold/Yellow	Orange/Red	Fructose	≤ ≩ 00
4H	p-n-p-β-D-glucoside	BGL	Yellow	Colorless	p-n-p-β-D-glucoside	≤¢0
2H	p-n-p-β-D-cellobioside	PCE	Yellow	Colorless	p-n-p-β-D-cellobioside	≤¢0
1H	Proline & Leucine-p- nitroanilide	PLN	Yellow	Colorless	Proline & Leucine-p- nitroanilide	≤¢0
41	p-n-p-phosphate	PHO	Yellow	Colorless	p-n-p-phosphate	ş 10
21	p-n-p-α-D-maltoside	PAM	Yellow	Colorless	p-n-p-α-D-maltoside	≤ ¢ 0
11	ONPG & p-n-p-α-D-galactoside	PGO	Yellow	Colorless	ONPG & p-n-p-α-D-galactoside	≤¢0
4J	Urea	URE	Aqua/Blue	Yellow/Green	Urea	≤ 50
2J	Esculin	ESC	Brown/Maroon	Clear/Tan	Esculin	≤ ¢5
1J	Arginine	ARG	Purple	Yellow/Gray	Arginine	≤ ≩ 00

9

Table 4 Quality (After 18	Table 4 Quality Control Chart for BBL Crystal tm GP ID System After 18–20 Hours Incubation from TSA II or Columbia Blood Agar	4 GP ID System A II or Columbia	Blood Agar	Table 5 Addition After 18	Table 5 Additional Quality Control Strains for BBL Crystal TM GP ID System After 18–20 Hours Incubation from TSA II or Columbia Blood Agar	BBL Crys A II or Co	tal™ GP ID Syste olumbia Blood Aç	m Jar		
Panel Location	Substrate	Code	Streptococcus pyogenes ATCC 19615	Panel Location	Substrate	Code	Staphylococcus epidermidis ATCC 12228	Bacillus brevis ATCC 8246	Enterococcus faecalis ATCC 19433	Staphylococcus xylosus ATCC 35033
4A	Fluorescent negative control	Ē		4A	Fluorescent negative control	FG	,	,	,	
ZA	4 MU-β-D-glucoside	FGC		74	4MtJ-8-D-alucoside	EGC		+	+	
14	L-valine-AMC	FVA	+					. :		
4B	L-phenylalanine-AMC	FPH	+	4		EV1		>		
28	4MU-\arapsilon-co-b-glucoside	FGS	+	4B	L-phenylalanine-AMC	FPH		+	+	
18	L-pyroglutamic acid -AMC	FΡΥ	+	2B	4MU-α-D-glucoside	FGS	*'	+	+	
4C	L-tryptophan-AMC	FTR	+	1B	L-pyroglutamic acid-AMC	FΡΥ		+	+	>
ž	L-arginine-AMC	FAR	+	4C	L-tryptophan-AMC	FTR		+	+	>
10	4MU-N-acetyl-β-D-glucosaminide	FGA		<u></u>			>	+		,
4	4MU-phosphate	FHO	+	ייי			>	+		
ą	4MU-β-D-glucuronide	FGN		5	4MU-N-acetyl-β-D-	FGA		+	+	
₽	L-isoleucine-AMC	FIS	+	!	giacosanimide					
4	Trehalose	TRE	+	4D	4MU-phosphate	EHO	+	>	>	+
2E	Lactose	LAC	+	2D	4MU-β-D-glucuronide	FGN				+
1	Methyl-α & β-glucoside	MAB	+	1D	L-isoleucine-AMC	FIS	-	^	-	
4	Sucrose	suc	+	4E	Trehalose	TRE			+	+
2F	Mannitol	MNT		2E	Lactose	LAC	+		+	+
Ŧ	Maltotriose	μt	+	;	Mothid ~ 9.8 alucorido					
4G	Arabinose	ARA		<u> </u>	ivietriyi-a & p-giucosiae				+	+
20	Glycerol	GLR	+	4F	Sucrose	SUC	+		+	+
5	Fructose	FRU	+	2F	Mannitol	MNT			+	+
4H	p-n-p-β-D-glucoside	BGL	>	1F	Maltotriose	TTM	+		+	*'
2H	p-n-p-β-D-cellobioside	PCE		4G	Arabinose	ARA				>
Ŧ	Proline & Leucine-p-nitroanilide	PLN	+	2G	Glycerol	GLR	+		+	+
4	p-n-p-phosphate	РНО	>	5	Enictoso	E	-	,	+	-
21	p-n-p-œ-D-maltoside	PAM	*'				÷		÷	÷
=	ONPG &	PGO		4H	p-n-p-b-U-glucoside	BGL		>	+	+
	p-n-p-α D-galactoside			2H	p-n-p-β-D-cellobioside	PCE			+	
4)	Urea	URE		Ħ	Proline & Leucine-p-nitroanilide	e PLN	>	>		
ה	Esculin	ESC		4	p-n-p-phosphate	ОНО	>	>	>	+
7	Arginine	ARG	>	21	p-n-p-α-D-maltoside	PAM	*'	>	+	*'
* = varia	= variable when tested from Columbia Blood Agar	Agar		=	ONPG & p-n-p-0:-D-galactoside	PGO	>	,	ı	>
				:	-					

10

¢

 2J
 Esculin
 ESC

 1J
 Arginine
 ARG

 * = variable when tested from Columbia Blood Agar

 \oplus

> +

> > + +

> > > +

> + > .

Urea

Ą

URE ESC ARG

> ï

8809911(0404).qxd 07/19/2004 3:18 PM Page 10

Table 6

Trials
Clinical
ystem
5
<u>-</u>
Σ
rystal™
BBLO
<u>ء</u> .
ntered
Encou
Ę
Frequent
ţ
ž
oecies
ŝ
ē
eactions
a B
Expecter

Organism	FCT	FCT FGC		FPH	FGS	FΡΥ	FTR	FAR	FGA	FHO	FGN	FIS T	RE L	AC N	AAB SL	IM DL	M TN	TT AI	FVA FPH FGS FPY FTR FAR FGA FHO FGN FIS TRE LAC MAB SUC MNT MTT ARA GLR FRU BGL PCE PLN PHO PAM PGO URE ESC ARO	FRU	BGL	PCE	PLN	РНО	PAM	I PGO	URE I	ESC	ARO
E. faecalis	ı	+	Т	+	+	+ + +	+	Т	+	Т	Т		+	>	+	+		>	+ + + + / + + (+) + / + -	+	+	+	I	>	+	(-) -	Ĵ	+	+
E. faecium	Т	+	I	+	(-)	+	+	Т	+	Т	Т	1	+	+	(+) / +		+	+	>	+	+	+	Т	(-)	+	(+)	I	+	+
S. aureus	I	I	I	I	I	(-)	>	+	I	+	I	-	(+	(+) / (+)		++	++		+	+	(+)	I	I	+	>	I	I	I	(+)
S. capitis	Т	I	I	I	Т	I	I	>	Т	Ĵ	I			I	- I - I	>	-		+	ŧ	T	I	T	>	I	I	I	I	(+
S. epidermidis	Т	Т	Т	Т	Т	Т	I	(-)	Т	+	Т	1		+	+		+		+	+	(-)	T	>	(+)	>	(-)	+	I	>
S. haemolyticus	I	I	I	I	I	+	I	(-)	>	- (-) - (-)	(-)	-	(+	>	+ (-) / (+)	>	+		+	(+	>	I	T	>	>)	I	I	>
S. hominis	Т	I	I	Т	Т	I	I	>	I	I	Т	-	>	>	>		+ +	- (+)	+	+)	I	>	>	>	Ĵ	+	I	I
S. agalactiae	I	I	>	+	(+	I	+	+	I	+	>	+ >		>	+	>	+		>	+	(+)	T	>	+	I)	I	I	>
S. bovis	Т	+	+	+	+	I	+	>	>	+ + (+) + ^ - ^ +	>	+	Ŧ	+	+	>	+		I	+	+	+	(+) + +	>	I	+	I	+	T
5. pneumoniae	I	>	+	+	+	I	+	+	+	I	Т	+	1	I	-	- I 			Ι	I	>	>	(+)	I	>	+	I	I	$\widehat{}$
S. pyogenes	Т	Т	+	+	>	+	+	+	Т	+ + (-) + -	-	+	+	+	+		+	т	(+)	+	- _	I	÷	+ (+)	Ĵ	I	I	>	>
KEY: + = 00% positive; (+)	ositiv	e; (+)		-89%	pos	itive;	V = 2	6-74	od %	sitive	-	11-	25%	positi	= 75–89% positive; V = 26–74% positive; (-) = 11–25% positive; - = ¢10% positive.	= ç10%	6 pos	itive.											

=00% positive; (+) = 75-89% positive; V = 26-74% positive; (-) = 11-25% positive; -= g10% positive. +

11

Table 7

Accuracy of Identification for Species Most Frequently Encountered in BBLCrystal™ GP ID System Clinical Trial

Organism	Number Tested	BBL Crystal Correct ID	BBL Crystal Correct W/Supplemental Tests	Total Correct
Corynebacterium species	29	29	0	29
Enterococcus cassaliflavus/gallinarum	14	0	14 ¹	14
Enterococcus faecalis	78	78	0	78
Enterococcus faecium	33	30	3	33
Micrococcus species	10	10	0	10
Staphylococcus aureus	88	85	3	88
Staphylococcus capitis	13	13	0	13
Staphylococcus epidermidis	87	87	0	87
Staphylococcus haemolyticus	23	23	0	23
Staphylococcus hominis	17	10	1	11
Streptococcus agalactiae	54	49	2	51
Streptococcus bovis	10	8	1	9
Streptococcus mitteri group	20	18	2	20
Streptococcus mitis group ²	23	8	1	9
Streptococcus pneumoniae	54	45	8	53
Streptococcus pyogenes	50	49	0	49
Other *	132	81	10	91
Grand Total	735	623	45	668

Key:

* = This category comprises all isolates where less than 10 were encountered in clinical trials.
 1 = Colony pigmentation is the sole supplemental test required to obtain correct identification.
 2 = As follow-up to this group's accuracy results, remedial actions were subsequently implemented to improve performance.

•••

Manufacturer / Výrobce / Producent / Fabrikant / Tootja / Valmistaja / Fabricant / Hersteller / Κατασκευαστής / Gyártó / Ditta produttrice / Gamintojas / Producent / Fabricante / Výrobca / Tillverkare

Use by / Spotřebujte do / Anvendes før / Houdbaar tot / Kasutada enne / Viimeinkäyttöpäivä / A utiliser avant / Verwendbar bis / Huepoµnía Aňčnç / Felhasználhatóság dátuma / Usare entro / Naudokite iki / Brukes før / Stosować do / Utilizar em / Použite do / Usar antes de / Använd före / YYYY-MM-DD / YYYY-MM (MM = end of month) / RRRR-MM-DD / RRR-MM (MM = konce méšice) ÁÅÅÅ-MM-DD / ÅÅÅÅ-MM (MM = slutning af måned) / JJJJ-MM-DD / JJJJ-MM (MM = einde maand) AAAA-KK-RP / AAAA-KK (KK = kuu löpp) VVVV-KK-PP / VVVV-KK (kuukauden loppuun mennessä) AAAA-KK-PP / AAA-KK (KK = kuu löpp) VVVV-KK-PP / VVVV-KK (kuukauden loppuun mennessä) AAAA-MM-JJ / AAAA-KM (MM = fin du mois) / JJJJ-MM-TT / JJJ-MM (MM = Monatsende) / EEEE-MM-HH / EEEE-MM (MM = téλoç rou µńya) / ÉÉÉÉ-HH-NN / ÉÉÉE-HH (HH = hónap utolsó napja) AAAA-MM-DD / ÅÅÅA-MM (MM = fine mese) / MMMM-MM-DD / ÅÅÅA-MM (MM = slutten av måneden) RRR-MM-DD / ÅÅÅA-MM (MM = fin do més) / RRR-MM-DD / RRR-MM (MM = fin del mes) / RRR-MM-DD / RRR-MM (MM = fin del mes) / AAAA-MM-DD / RRR-MM (MM = fin del mes) / AAAA-MM-DD / RRR-MM (MM = fin del mes) / AAAA-MM-DD / RRR-MM (MM = fin del mes) / AAAA-MM-DD / RRR-MM (MM = sluttet på månaden)

REF

Catalog number / Katalogové číslo / Katalognummer / Catalogusnummer / Kataloogi number / Tuotenumero / Numéro catalogue / Bestellnummer / Αριθμός καταλόγου / Katalógusszám / Numero di catalogo / Katalogo numeris / Numer katalogowy / Número do catálogo / Katalógové číslo / Número de catálogo

EC REP

Authorized Representative in the European Community / Autorizovaný zástupce pro Evropskou unii / Autoriseret repræsentant i EU / Erkend vertegenwoordiger in de Europese Unie / Volitatud esindaja Euroopa Nõukogus / Valtuutettu edustaja Euroopan yhteisössä / Représentant agréé pour la C.E.E. / Autorisierte EG-Vertretung / Eξουσιοδοτημένος αντιπρόσωπος στην Euρωπαϊκή Κοινότητα / Hivatalos képviselet az Európai Unióban / Rappresentante autorizzato nella Comunità europea / Įgaliotasis atstovas Europos Bendrijoje / Autorisert representant i EU / Autoryzowane przedstawicielstwo w Unii Europejskiej / Representante autorizado na União Europeia / Autorizovaný zástupca v Európskom spoločenstve / Representante autorizado en la Comunidad Europea / Auktoriserad representant i EU

IVD

In Vitro Diagnostic Medical Device / Lékařské zařízení určené pro diagnostiku in vitro / In vitro diagnostik medicinsk anordning / Medisch hulpmiddel voor in vitro diagnoste / In vitro diagnostik meditsiiniaparatuur / Lääkinnällinen in vitro -diagnostiikkalaite / Dispositif médical de diagnostic in vitro / Medizinisches In-vitro-Diagnostikum / In vitro διαγγωστική ατοική συσκευή / In vitro diagnostika i orvosi eszköz / Dispositivo medico diagnostic in vitro / In vitro diagnostika j In vitro diagnostika medisinsk utstyr / Urządzenie medyczne do diagnostyki in vitro / Dispositivo médico para diagnóstico in vitro / Medicínska pomôcka na diagnostiku in vitro / Dispositivo médico de diagnostico in vitro / Medicínska nordning för in vitro-diagnostik

Temperature limitation / Teplotní omezení / Temperaturbegrænsning / Temperatuurlimiet / Temperatuuri piirang / Lämpötilarajoitus / Température limite / Zulässiger Temperaturenbereich / Όριο θερμοκρασίας / Hőmérsékleti határ / Temperatura limite / Laikymo temperatūra / Temperaturbegrensning / Ograniczenie temperatury / Limitação da temperatura / Ohraničenie teploty / Limitación de temperatura / Temperaturbegränsning

LOT

Batch Code (Lot) / Kód (číslo) šarže / Batch kode (Lot) / Chargenummer (lot) / Partii kood / Erăkoodi (LOT) / Code de lot (Lot) / Chargencode (Chargenbezeichnung) / Κωδικός παρτίδας (Παρτίδα) / Tétel száma (Lot) / Codice del lotto (partita) / Partijos numeris (Lot) / Batch-kode (Serie) / Kod partii (seria) / Código do lote (Lote) / Kód série (šarža) / Código de lote (Lote) / Satskod (parti)

i

Consult Instructions for Use / Prostudujte pokyny k použití / Læs brugsanvisningen / Raadpleeg gebruiksaanwijzing / Lugeda kasutusjuhendit / Tarkista käyttöohjeista / Consulter la notice d'emploi / Gebrauchsanweisung beachten / Συμβουλευτείτε τις οδηγίες χρήσης / Olvassa el a használati utasítást / Consultare le istruzioni per l'uso / Skaitykite naudojimo instrukcijas / Se i bruksanvisningen / Zobacz instrukcja użytkowania / Consulte as instruções de utilização / Pozri Pokyny na používanie / Consultar las instrucciones de uso / Se bruksanvisningen

BD

Becton, Dickinson and Company 7 Loveton Circle Sparks, Maryland 21152 USA 800-638-8663

ATCC is a trademark of the American Type Culture Collection. BD, BD Logo, BBL, BBL Crystal, and Trypticase are trademarks of Becton, Dickinson and Company. © 2004 BD.